



**PRODUCT CODE**

NC10030

*Access Point Medical™ Deluxe Nebulizer*

**Product Overview** The Access Point Medical™ Deluxe Nebulizer provides a source of compressed air for medical purposes for use in home health care. The Deluxe Nebulizer is used in conjunction with a pneumatic nebulizer to produce a fine aerosol mist of medication for respiratory therapy, for both children and adults suffering from respiratory disorders.

- Features**
- Compact and simple to use
  - Attractive retail packaging for easy display
  - Includes MicroMist™ disposable nebulizer kit

**Warranty** 5-year limited warranty

|  |                            |
|--|----------------------------|
| <b>ELECTRIC REQUIREMENTS</b>           | AC 115V/60Hz               |
| <b>COMPRESSOR<br/>FREE AIR FLOW</b>    | ≥15L/min                   |
| <b>SOUND LEVEL</b>                     | ≤ 64dB (A)                 |
| <b>WEIGHT</b>                          | 3.3 lbs                    |
| <b>DIMENSIONS</b>                      | 9.5" (L) 4.8"(W) × 6.8"(H) |
| <b>ELECTRIC SAFETY</b>                 | UL60601-1                  |
| <b>WORKING COMPRESSOR<br/>PRESSURE</b> | > 10 psi, Max ≥ 30 psi     |
| <b>MAXIMUM RUNNING TIME</b>            | 30 Minutes                 |
| <b>NORMAL AMBIENT RANGES</b>           | 50° F ~ 104° F             |
| <b>UNITS PER CASE</b>                  | 6                          |

MicroMist™ is a trademark or registered trademark of Hudson RCI.